



# Biohazard Waste Cans

To store contaminated laundry as well as other regulated waste



# Biohazard Waste Cans

## Waste Disposal Safety Cans - Biohazard Waste Cans

Justrite Biohazard Waste Containers meet OSHA 29 CFR 1910.1030 for the segregation of infectious waste. Use them to store contaminated laundry as well as other regulated waste; however, not for the storage of sharps.



- Safe, temporary collection of biological waste
- Meets OSHA 29 CFR 1910.1030 for segregation of infectious waste
- Foot-operated, self-closing lid
- Lessen disruptive closure noise with SoundGard™ models
- Approved by FM Approvals and TUV

- For the safe and temporary collection of biological, medical, contaminated or hazardous waste and laundry in industries like medical, laboratory and pharma.

Recommended use



# Waste Disposal Safety Cans - Biohazard Waste Cans

## Biohazard Waste Cans

### General

Brand	Justrite
Category	Biohazard Waste Cans
Product Number	0593
Product Family	Safety Cans & Containers
Product Group	Waste Disposal Safety Cans
Industry	Education, Health Care, Laboratories, Industrial, Manufacturing, Facility Management
Product Name	Biohazard Waste Cans
Recommended use	For the safe and temporary collection of biological, medical, contaminated or hazardous waste and laundry in industries like medical, laboratory and pharma.
Shipping Location	Achterzeedijk 57, Unit 15, Barendrecht 2992SB, The Netherlands
Update trigger	20250521_0815-ecom-fp

### Technical Specifications

Material	Steel
Material Group	Steel

## Biohazard Waste Can - White - Foot-Operated - Self-Closing - SoundGard

Art. nr.

JCN05935WH

## Biohazard Waste Can - White - Foot-Operated - Self-Closing

Art. nr.

JCN05910WH

JCN05930WH

## Biohazard Waste Can - Red - Foot-Operated - Self-Closing

Art. nr.

JCN05930RD

JCN05910RD

## Accessories

5920

### **Biohazardous Waste Bags**

For the safe,  
temporary collection  
of biological waste.  
Red colour is easily  
identifiable as  
biohazard bags.  
14 microns thick,  
available in packs of  
100.